

Medicraft **Hill-Rom®**

SPARE PART REFERENCE GUIDE

FET5000



Revision 1


Written and Constructed by: *Mark Griffith*
Authorised by: *Robert Crawford*

Last updated: 28 April 2009

Section 1 : Common Electrical Parts Linak



Description	Part Number	Supersession	Picture
Main Lift Actuator	710460	Med-SP-5 96	
Cables to suit - 710460	710466		
Backrest Actuator (New)	710461	712021	
QR cable suit – 710461	710503		
Mini Fit suit - 710461	710465		
Backrest Actuator (Old) <i>Note: Discontinued</i> <i>Refer to: 710461</i>	712021		N/A Discontinued Product
QR cable suit - 712021	710498	710503	

Knee Break Actuator (Hard Wired Pull Type)	710688	710459	A black, rectangular device with a coiled black cable, resting on a blue and white striped fabric with the text 'Hill-Rom' visible.
Knee Break Actuator (Mini Fit)	710459		A black, rectangular device with a coiled black cable, resting on a blue and white striped fabric with the text 'Hill-Rom' visible.
ACP Electric Tilt (Large)	710456	Med-SP-5 91	A white, rectangular device with a coiled white cable, resting on a blue and white striped fabric with the text 'Hill-Rom' visible.
Mini ACP Lockout - <i>optional</i>	710672	Med-SP-5 93	A white, rectangular device with a coiled white cable, resting on a blue and white striped fabric with the text 'Hill-Rom' visible.
Control Box CB12	710457	Med-SP-5 95	A white, rectangular device with a coiled white cable, resting on a blue and white striped fabric with the text 'Hill-Rom' visible.
Power Cable IEC - CB12	710463	Med-SP-5 179	A black cable with a red and black connector, resting on a blue and white striped fabric with the text 'Hill-Rom' visible.
Power Cable IEC - Wall	710464	Med-SP-5 143	A black cable with a red and black connector, resting on a blue and white striped fabric with the text 'Hill-Rom' visible.

Handset 5 channel HB75	710458	Med-SP-5 98	
Handset 4 channel HB74	710684	Med-SP-5 89	
Handset Extension Cable	710462	Med-SP-5 197	
Battery Backup set	711253	Med-SP-5 94	



Section 2: Common Electrical Parts Dewert

Description	Part Number	Supersession	Picture
Main Lift Actuator	710695		
Backrest Actuator (standard)	710696		
Quick Release Cable Suit 710696	710498		
Backrest Actuator (sliding)	710697		
Knee Break Actuator	710698		
ACP Electric Tilt (Large)	710700		

Control Box	710699		
Handset 5 Channel	710701		
Handset Extension Cable	710147		
Power Cable IEC – Wall	711052		
Power Cable IEC- CB ‘B’ type 2.5m	710705		
Battery Back Up Set	712699		

Section 3: Common Sleep Surface and Frame Parts

Description	Part Number	Supersession	Picture
ICO backrest section panel standard (665x770mm)	710178		 A white plastic backrest section panel with a grid of 12 circular holes arranged in 3 rows and 4 columns.
ICO fixed section panel standard (200x770mm)	710179		 A white plastic fixed section panel with a single row of 4 circular holes.
ICO middle section panel standard (365x770mm)	710180		 A white plastic middle section panel with a single row of 4 circular holes.
ICO foot section panel standard (565x770mm)	710181		 A white plastic foot section panel with a grid of 12 circular holes arranged in 3 rows and 4 columns.
Plastic ICO handles	710182		 A white plastic handle with a D-shaped cutout.
Sheet Metal Backrest Section Panel	712293		 A white sheet metal backrest section panel with a grid of 12 circular holes arranged in 3 rows and 4 columns.

Sheet Metal Fixed Section Panel	710581		
Sheet Metal Middle Section Panel	710582		
Sheet Metal Foot Section Panel	712294		
Swing Telescopic (Electric Tilt)	710516 x 1 710509 x 1 710664 x 1 710517 x 1	Med-SP-16 70	
Swing 'A' Frame	710516	Med SP 16 38	
Height Lever	710204	Med-SP-4 31	

Hydraulic Damper/ Gas Strut (Linak systems made before '08)	711246		
Horizontal Safety Siderail Set	710748		
Horizontal Safety RH Siderail	710750		
Horizontal Safety LH Siderail	710751		
Horizontal Safety Siderail Set (Slimline Design) Optional	713164		
Horizontal Safety Siderail RH (Slimline Design) Optional	713163		
Horizontal Safety Siderail LH (Slimline Design) Optional	713162		

Side Rail release assembly	710851	Med-SP-15 79	
Wall protection buffers	712124 712125	(Grey) (Cream)	


Section 4: Common Brake Assembly Parts

Description	Part Number	Supersession	Picture
150mm Twin Wheel Castor Set Optional	710607		
150mm Twin Wheel Directional Lock Castor Optional	710988		
150mm Twin Wheel Total Lock Castor Optional	710987		
150mm Single Wheel Castor Set	710604		
150mm Single Wheel Directional Lock Castor	710414	Med-SP-914	
150mm Single Wheel Total Lock Castor	710413		

200mm Single Wheel Directional Lock Castor Optional	710985		
200mm Single Wheel Total Lock Castor Optional	710986		
Wheel Base Hex Bar – Head	710197		
Wheel Base Hex Bar - Foot	710191		
Wheel Base Central Lock Link	710526		
Transmission Lever (Tente)	711520		


Brake Pedal Left Hand	710202		
Brake Pedal Right Hand	710199		

Linak Actuator Identification Figure 1.1




Note: Picture is from an under bed vantage

Dewert Actuator Identification Figure 2.1



Note: Picture is from an under bed vantage.

ICO Panel Identification Figure 3.1



Please Note: ICO stands for “Infection Control Option” and is a sleep surface panel especially designed to allow for ease of cleaning and to reduce the spread of infection. The panels are made with injection moulded plastic with a high percentage of non-recycled material.

Notations:

Part Number effected:	Notations:

Medicraft **Hill-Rom**

Customer Service & Technical Advice

Free Call: 1800 445 576

Ph: 02 9569 0255 Fax: 02 9569 8394